

**Co-Immunoprecipitation of germline proteins from adult *C. elegans***  
**August 2014 Christian's Method**

Part 1: Grow worms & make popcorn (Dustin's protocol w/ modifications)

1. Make worm popcorn from synchronized adults
  - 1.1. Day 1 Grow up worms to synchronize
    - 1.1.1. Add 25mL of HB101 bacteria to 1 liter of complete S media.
    - 1.1.2. Wash worms off of 40 recently starved plates and add to the liter of S media and bacteria.
    - 1.1.3. Divide S media + bacteria + worms into 4x 250mL glass flasks. Shake at 220rpm at 20C for 3 days.
  - 1.2. Day 4 late: Harvest embryos and hatch on unseeded plates
    - 1.2.1. Pour worms into 1 liter pear funnel and let settle.
    - 1.2.2. Drip into 50ml conical tube.
    - 1.2.3. Wash with water 3X
    - 1.2.4. Throw worms into bleach solution. (~3 minutes then rinse 3X with M9)
      - 1.2.4.1. 3 parts water, 1 part bleach, 0.1 part 10M NaOH (50ml – 36mL water, 12.6mL bleach, 1.4mL 10N NaOH)
    - 1.2.5. Let embryos hatch on about 20 unseeded (no food) plates.
  - 1.3. Day 5 mid-day: Feed L1's
    - 1.3.1. Add 40mL of HB101 bacteria to 1 liter of complete S media.
    - 1.3.2. Wash hatched L1s and add to the liter of S media and bacteria
    - 1.3.3. Divide S media + bacteria + worms into 4x 250mL glass flasks. Shake at 220rpm at 20C for 24 hours.
  - 1.4. Day 6 early: Reduce temp
    - 1.4.1 Reduce temp from 20C to 18C.
    - 1.4.2 Otherwise worms will be too old by harvesting time on Day 8.
  - 1.5. Day 8 early Harvest adults and make popcorn
    - 1.5.1. Inspect worms on dissection scope to make sure most are adults (around 68 hours after feeding L1s).
    - 1.5.2. Pour worms into 1 liter pear funnel and let settle
    - 1.5.3. Pre-chill 50ml conical tubes in dry ice
    - 1.5.4. Collect and wash 2X in water.
    - 1.5.5. Wash 1X in Buffer H (lysis buffer) with NO detergent!! (see part 3 for recipe)
    - 1.5.6. Pour liquid nitrogen into 50ml pre-chilled tube
    - 1.5.7. Put pelleted worms in syringe and drip worms into liquid nitrogen
    - 1.5.8. Put worm popcorn in 50ml conical tubes and freeze in -80C

This protocol usually nets somewhere around 10-20mL of packed popcorn depending on how healthy your strain is and the # of adults you have at bleaching stage. I've had terrible luck using worms from NGM plates and ALWAYS do liquid cultures now.

### Part 2: Make worm powder

- 2.1. Pulverize worms in Retsch Mixer Mill MM-400 under cryogenic conditions
  - 2.1.1. Chill 2x 25mL grinding jars w/ balls inside in liquid nitrogen until bubbling stops
  - 2.1.2. Transfer around ~15mL of popcorn to each tube, leaving room for ball to move around. DO NOT OVERFILL.
  - 2.1.3. NOTES: You must ALWAYS use both jars, similar to balancing centrifuge. Both jars MUST have popcorn in them, ball will damage an empty jar.
  - 2.1.4. After filled with popcorn, return ball to jar, seal, and immerse in liquid nitrogen.
  - 2.1.5. Keep jars in liquid nitrogen when not milling.
  - 2.1.6. Pulverize worms for 2min @ 25Hz, 3 times. Chill on liquid nitrogen between milling sessions until bubbling stops.
  - 2.1.7. Carefully remove powder from jars into small box or conical tube.
  - 2.1.8. Store at -80C until ready to perform IP (can be stored indefinitely).

This usually produces somewhere in the ~5g range of powder, enough for a few IPs.

### Part 3: Conjugate dynabeads to desired antibody (Bungo's protocol)

I always cross-link my antibodies to Invitrogen Dynabeads. This allows you to use the same antibody for Western blot without having to worry about heavy & light chains showing up. For a single Co-IP, I make 100uL of beads (50uL each for sample/negative control) w/ 25ug of antibody (12.5ug/IP).

Use Protein A Dynabeads for rabbit antibodies and Protein G for mouse antibodies.

Working IP antibodies:

- mouse GFP (Roche, 11814460001)
- mouse FLAG M2 (Sigma, F3165)
- rabbit A-14 Myc (Santa Cruz Biotech, SC-789)

- 3.1. Bind 25ug of antibody to 100uL of Dynabeads.
  - 3.1.1. Resuspend 100uL slurry of Dynabeads.
  - 3.1.2. Concentrate in magnet, remove sup.
  - 3.1.3. Wash twice in 1mL of 0.1M Na-phosphate buffer pH 7.0
    - 3.1.3.2. 500mL 0.1M Na phosphate pH 7.0 (58mL 0.5M Na<sub>2</sub>HPO<sub>4</sub>, 21mL 1M NaH<sub>2</sub>PO<sub>4</sub>, 421mL ddH<sub>2</sub>O)
  - 3.1.4. Add 25ug of antibodies and same volume of 0.1M Na phosphate buffer.
    - 3.1.4.2. e.g. if 25ug of FLAG M2 is 25uL, add 25uL of 0.1M Na phosphate
  - 3.1.5. Shake with gentle agitation on vortex machine for 25min.
- 3.2. Crosslink antibody to beads
  - 3.2.1. Concentrate beads, remove sup.
  - 3.2.2 Wash twice in 1ml of 0.1M Na-phosphate buffer pH 7.0 with 0.01% Tween-20 (Tween helps with protein stability in storage of Dynabeads)
  - 3.2.3. Wash twice with 1ml of 0.2M triethanolamine, pH 8.2 (Dynabeads now have immobilized antibody)
  - 3.2.4. Resuspend Dynabeads in 1ml of 20mM DMP (in common -20°C freezer, Dimethyl Pimelimidate, D8388, Sigma) in 0.2M triethanolamine, pH 8.2 (note that this crosslinking solution must be prepared immediately before use)
    - 3.2.4.2. Use 5.18mg DMP in 1mL 0.2M triethanolamine
  - 3.2.5. Incubate on rotation mixer for 30 minutes at RT (DO NOT OVER CROSSLINK)
- 3.3. Stop the crosslinking reaction
  - 3.3.1. Concentrate beads, remove sup.
  - 3.3.2. Resuspend beads in 1mL of 50mM Tris-HCl pH=7.5
  - 3.3.3. Incubate with rotation mixing for 15min at RT.
  - 3.3.4. Wash beads 3X in 1mL PBST (0.1% Tween-20)
  - 3.3.4 Store beads in original volume but w/ PBST (e.g. 100uL PBST) at 4C. Do not freeze.
  - 3.3.5 Beads can be stored for around ~1 month at 4C without any reduction in efficiency. No need to equilibrate prior to use.

#### Part 4: Lyse Worms & perform IP

Lysis buffer (Buffer H): 50mM Hepes pH 8.0, 150mM KCl, 2mM MgCl<sub>2</sub>, 15% Glycerol, 0.1mM EDTA pH 8.0, 0.5mM EGTA-KOH pH 8.0, 0.1% NP-40

Protease inhibitors (working concentrations): 0.1 mM AEBSF, 5 mM benzamidine, 1:200 aprotinin, Roche® Complete Mini tablets w/o EGTA

Phosphatase inhibitors (working concentrations): 1mM Na<sub>4</sub>P<sub>2</sub>O<sub>7</sub>, 2mM Na-β-glycerophosphate, 0.1mM Na<sub>3</sub>VO<sub>4</sub>, 5mM NaF

2X Sample Buffer (2XSB): 4% SDS, 20% Glycerol, 120mM Tris-HCl pH 6.8, 0.02% bromophenol blue (add 2-mercaptoethanol [B-Me] to 5% prior to use)

#### 4.1. Lyse worms

- 4.1.1. Weigh 1g of worm powder for each strain you plan to use into 15mL conical tube (always do untagged negative control).
- 4.1.2. Resuspend in 5mL ice cold Buffer H w/ protease + phosphatase inhibitors.
- 4.1.3. Rotate 3min @ 4C or until re-suspended well.
- 4.1.4. Sonicate 2x30sec at medium amplitude with 1min rest on ice.
- 4.1.5. Spin @ ~50,000xg (max speed) @4C for 20' in JA-20 rotor.
- 4.1.5.2. Can do optional ~15min ultracentrifuge spin @ 4C. This helps get rid of lipids and seems to help with IP, but not necessary.
- 4.1.6. Take sup, move to new 15mL conical avoiding top lipid layer.

#### 4.2. Perform IP

- 4.2.1. Take 50uL 'input' sample for each strain, add 50uL 2X SB+B-Me, boil 3 min, store at -20C.
- 4.2.2. Add 50uL of Dynabeads with conjugated antibody.
- 4.2.2.2. Note this is usually not enough to immuno-deplete, could increase to get more protein.
- 4.2.3. Rotate at 4C for 2-3hours.
- 4.2.4. Concentrate beads in magnet.
- 4.2.5. Take 50uL of sup as 'unbound' sample for each strain, add 50uL 2X SB+B-Me, boil 3 min, store at -20C.
- 4.2.6. Remove sup, wash beads 4X with 1mL ice cold Buffer H.
- 4.2.7. Remove as much sup as possible.
- 4.2.8. Boil in 50uL 2X SB+B-Me to elute IP'd proteins.
- 4.2.9. Concentrate in magnet, take sup (eluate) to new tube (removing beads).

### Part 5: Analyze IP'd proteins by Western blot

See Western protocol.

I usually run 10% of IP sample (5uL) to analyze for the presence of IP'd protein (e.g. protein with GFP/FLAG tag)

To analyze for presence of Co-IP'd proteins, I run significantly more, up to 50% (25uL) of IP sample, depending on stoichiometry of complex.

For inputs/unbound samples, I run 10uL each for Western.